


Be Extra Smart,
Be Extra Refined

X8 Patient Monitor



X8
Patient Monitor



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com


© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

Build Your Confidence on Quality Care

Based on its expandability of monitoring parameters, and stunning ultra-slim design, the X8 fulfills your primary clinical requirements in various scenarios, including Emergency Rooms, General Wards, Rehabilitation Department, Cardiac Units, and In-hospital Transfer, covering all patient types, realizing from neonatal to adult monitoring.



8"




No-fan Design



Pitch Tone



Dual Alarm Lights



Shortcut Menu



Wi-Fi



USB Disk



Operating Time ≥8h



Barcode Scanner

Standard Parameters:

3/5-lead ECG, RESP, EDAN SpO₂, EDAN NIBP, TEMP


Optional Parameters:

6/12-lead ECG, EDAN G2 CO₂

Proprietary Algorithms & Technologies

EDAN G2 CO₂ (Sidestream)

- Superior water trap design for accurate monitoring.
- iCARB™ algorithm with intelligent CO₂ pseudo wave identification technology.
- Sampling rate as low as 50ml/min.
- Accessories for all patient types.



ECG

- Customizable 6-lead placement for more ECG waves.
- Automatic lead type detection.
- Industry leading iSEAP™ algorithm with auto-detection of 33 types of arrhythmias.
- SEMIP™ algorithm with 208 ECG findings over age/gender diversities.



NIBP

- Dual dust filter design makes no blockage inside and provides accurate NIBP readings.
- Unique cleaning mode for routine maintenance.
- iCUFS™ algorithm with smart deflation technology.

SpO₂

- iMAT™ algorithm with outstanding motion resistance and low perfusion resistance performance.
- Reference reading of Perfusion Index (PI) from 0 to 10 according to perfusion changes.
- Simultaneous measurements of SpO₂ and NIBP of the same limb.



A single piece of patient data contains:

1200

NIBP

240 h

Trend Review

200

Alarm Review

48 h

Waves